



Protection 24/7

SaveStations can now be placed outdoors in high traffic public places, providing the public with unprecedented access to life-saving AEDs.

REMOTE MONITORING

Sensor monitoring the presence and ready status of the AED and the cabinet climate.

CONNECTIVITY

Wireless communication.

REAL-TIME REPORTING

Immediate notifications via email/SMS.

CABINET LIGHTING

Integrated light sensor to illuminate the cabinet.

PHOTO CONTROL

Infra-red photo control when cabinet is opened.

LCD SCREEN

Operational status visual at all times for manual checks via integrated LCD screen.

POWER

24V DC convertor - SELV (Safety Extra-Low Voltage).

ALARM

Audible alarm and visual warnings triggered on door opening (70 dB at 1 metre, and flashing red LEDs).

TEMPERED GLASS

Impact-resistant transparent polycarbonate door with surface treatment to minimize greenhouse effect.

REGULATED HEATING AND VENTILATION

Maintains positive temperature of AED with ventilator fan and convector heating device.

BACK-UP POWER

Rechargeable back-up battery for essential functions in the event of a power supply failure.



info@actionfirstaid.ca / 1-866-347-7824

Action First Aid.ca
Empowering People to **TAKE ACTION**